

## Pseudocowpox Virus (VR634) PCPV\_gp006 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX560

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Pseudocowpox Virus (VR634) PCPV\_gp006 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                    |
|----------------|--------------------|
| Target Protein | PCPV_gp006         |
| Species        | Pseudocowpox Virus |
| Strain         | VR634              |
| Protein Class  | Virology           |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                   |
|-------------------|-----------------------------------|
| Target Protein    | PCPV_gp006                        |
| Full Name         | Ankyrin/F-box protein             |
| Alternative Names | PCPV_gp006; Ankyrin/F-box protein |
| Gene ID           | <a href="#">8797207</a>           |
| UniProt ID        | <a href="#">D3IZF1</a>            |