

Volepox virus (CA) BI079_gp140 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX57

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Volepox virus (CA) BI079_gp140 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	BI079_gp140
Species	Volepox virus
Strain	CA
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	BI079_gp140
Full Name	DNA polymerase processivity factor
Alternative Names	BI079_gp140; DNA polymerase processivity factor; VPXV-CA-140
Gene ID	29064087
UniProt ID	A0A1C9KCH0