

Parapoxvirus of Red Deer (HL953) SB87_gp008 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX687

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Parapoxvirus of Red Deer (HL953) SB87_gp008 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SB87_gp008
Species	Parapoxvirus of Red Deer
Strain	HL953
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SB87_gp008
Full Name	putative EEV phospholipase
Alternative Names	SB87_gp008; putative EEV phospholipase
Gene ID	22647408
UniProt ID	A0A0A7MEL7