

## Volepox virus (CA) BI079\_gp127 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX70

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Volepox virus (CA) BI079\_gp127 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

Target Protein	BI079_gp127
Species	Volepox virus
Strain	CA
Protein Class	Virology

### Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### Target

Target Protein	BI079_gp127
Full Name	viral membrane associated protein
Alternative Names	BI079_gp127; viral membrane associated protein; VPXV-CA-127
Gene ID	<a href="#">29064074</a>
UniProt ID	<a href="#">A0A1C9KCD0</a>