

Parapoxvirus of Red Deer (HL953) SB87_gp032 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX717

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Parapoxvirus of Red Deer (HL953) SB87_gp032 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SB87_gp032
Species	Parapoxvirus of Red Deer
Strain	HL953
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SB87_gp032
Full Name	hypothetical protein
Alternative Names	SB87_gp032; hypothetical protein
Gene ID	22647432
UniProt ID	A0A0A7M9Q2