

Squirrel Parapoxvirus (Red squirrel UK) A28L Recombinant Protein

Made to Order

Cat: MPYF-1022-KX759

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Squirrel Parapoxvirus (Red squirrel UK) A28L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	A28L
Species	Squirrel Parapoxvirus
Strain	Red squirrel UK
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	A28L
Full Name	IMV surface protein
Alternative Names	A28L; SQPV_1170; IMV surface protein
Gene ID	18158326
UniProt ID	Q1HTR8