

Squirrel Parapoxvirus (Red squirrel UK) J1R Recombinant Protein

Made to Order

Cat: MPYF-1022-KX781

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Squirrel Parapoxvirus (Red squirrel UK) J1R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	J1R
Species	Squirrel Parapoxvirus
Strain	Red squirrel UK
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	J1R
Full Name	subunit of the poxvirus multiprotein entry-fusion complex
Alternative Names	J1R; SQPV_0590; subunit of the poxvirus multiprotein entry-fusion complex
Gene ID	18158409
UniProt ID	U3UB99