

Squirrel Parapoxvirus (Red squirrel UK) G1L Recombinant Protein

Made to Order

Cat: MPYF-1022-KX793

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Squirrel Parapoxvirus (Red squirrel UK) G1L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	G1L
Species	Squirrel Parapoxvirus
Strain	Red squirrel UK
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	G1L
Full Name	metalloproteinase G1-MEROPS peptidase family M44 (clan ME)
Alternative Names	G1L; SQPV_0440; metalloproteinase G1-MEROPS peptidase family M44 (clan ME)
Gene ID	18158396
UniProt ID	U3UB96