

Squirrel Parapoxvirus (Red squirrel UK) D2L Recombinant Protein

Made to Order

Cat: MPYF-1022-KX815

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Squirrel Parapoxvirus (Red squirrel UK) D2L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	D2L
Species	Squirrel Parapoxvirus
Strain	Red squirrel UK
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	D2L
Full Name	virion core protein / unknown function
Alternative Names	D2L; SQPV_0730; virion core protein / unknown function
Gene ID	18158352
UniProt ID	U3UB19