

Squirrel Parapoxvirus (Red squirrel UK) SQPV_0350 Recombinant Protein Made to Order

Cat: MPYF-1022-KX866

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Squirrel Parapoxvirus (Red squirrel UK) SQPV_0350 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SQPV_0350
Species	Squirrel Parapoxvirus
Strain	Red squirrel UK
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SQPV_0350
Full Name	conserved hypothetical pox protein
Alternative Names	SQPV_0350; conserved hypothetical pox protein
Gene ID	18158361
UniProt ID	U3UBC4