

Anti-Skunkpox virus SKPV-WA-160 Monoclonal Antibody (MP-K853) (Mouse IgG)

Cat: MPYF-1022-KX919

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Skunkpox virus SKPV-WA-160. The mAb with highest affinity will be selected. The selected antibody recognizes Skunkpox virus SKPV-WA-160. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	SKPV-WA-160
Species Reactivity	Skunkpox virus
Strain	WA
Cross Reactivity	SKPV-WA-160
Specificity	This antibody recognizes Skunkpox virus SKPV-WA-160.
Antibody Isotype	Mouse IgG
Clone	MP-K853
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	SKPV-WA-160
Gene Name	SKPV-WA-160
Introduction	Skunkpox viruses are dsDNA viruses that reproduce in the cytoplasm of infected eukaryotic cells. Colfax, Washington (USA), is where the SKPV-WA was isolated in 1978. It is unknown, though, if these infections were linked to any serious illnesses. In support of prior genome hybridization investigations, early DNA sequencing of individual genes revealed that these NAOVs were more divergent than the Old World Afri

can and Eurasian orthopoxviruses and that, despite appearing to constitute a different clade, they were likely still orthopoxviruses.

Alternative Names	BIZ96_gp160; semaphorin; SKPV-WA-160
Official Symbol	Semaphorin
Gene ID	29057738
UniProt ID	A0A1C9KBV9