

Anti-Bovine Papular Stomatitis Virus BPSVgORF045 Monoclonal Antibody (MP-K1891) (Mouse IgG)

Cat: MPYF-1122-KX1091

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Bovine Papular Stomatitis Virus BPSVgORF045. The mAb with highest affinity will be selected. The selected antibody recognizes Bovine Papular Stomatitis Virus BPSVgORF045. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	BPSVgORF045
Species Reactivity	Bovine Papular Stomatitis Virus
Strain	BV-AR02
Cross Reactivity	BPSVgORF045
Specificity	This antibody recognizes monkeypox virus (MPXV) BPSVgORF045.
Antibody Isotype	Mouse IgG
Clone	MP-K1891
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	BPSVgORF045
Gene Name	BPSVgORF045
Introduction	Similar to other poxviruses, the BPSV genome contains a crucial central coding region delimited by two inverted terminal repeat (ITR) regions.

Alternative Names	BPSVgORF045; late transcription factor VLTF-1; viral transcriptional activator of late genes; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Official Symbol	late transcription factor VLTF-1
Gene ID	2947939
UniProt ID	Q6TVE3