

Anti-Orf Virus (Sore Mouth Infection) ORFVgORF100 Monoclonal Antibody (MP-K2026) (Mouse IgG)

Cat: MPYF-1122-KX1226

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Orf Virus (Sore Mouth Infection) ORFVgORF100. The mAb with highest affinity will be selected. The selected antibody recognizes Orf Virus (Sore Mouth Infection) ORFVgORF100. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	ORFVgORF100
Species Reactivity	Orf Virus (Sore Mouth Infection)
Strain	OV-SA00
Cross Reactivity	ORFVgORF100
Specificity	This antibody recognizes monkeypox virus (MPXV) ORFVgORF100.
Antibody Isotype	Mouse IgG
Clone	MP-K2026
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	ORFVgORF100
Gene Name	ORFVgORF100
Introduction	ORFV, a typical member of the Poxviridae family, is a large oval double-stranded DNA virus. The ORFV genome is approximately 140 kb and contains 132 genes. ORF

Vs include highly variable terminal regions and relatively conserved central regions. Terminal regions are involved in virulence, pathogenesis, and immune evasion.

Alternative Names	ORFVgORF100; ORF100 intermediate transcription factor VITF-3
Official Symbol	ORF100 intermediate transcription factor VITF-3
Gene ID	2947715
UniProt ID	Q6TVM0