

---

## Anti-Pseudocowpox Virus PCPV\_gp124 Monoclonal Antibody (MP-K2138) (Mouse IgG)

**Cat: MPYF-1122-KX1338**

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Pseudocowpox Virus PCPV\_gp124. The mAb with highest affinity will be selected. The selected antibody recognizes Pseudocowpox Virus PCPV\_gp124. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

### Product Description

Target	PCPV_gp124
Species Reactivity	Pseudocowpox Virus
Strain	VR634
Cross Reactivity	PCPV_gp124
Specificity	This antibody recognizes monkeypox virus (MPXV) PCPV_gp124.
Antibody Isotype	Mouse IgG
Clone	MP-K2138
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

### Target

Target	PCPV_gp124
Gene Name	PCPV_gp124
Introduction	PCPV belongs to the family of poxviruses, which includes several members of the virus with economic impact. PCPV has a linear double-stranded DNA genome of 145 kb in length with high GC content and 134 open reading frames (ORFs). The PCPV linear genome contains central and terminal regions. The central region is

responsible for expressing structural proteins and factors related to viral proliferation while the two terminal regions are closely related to genes related to virulence factors.

---

Alternative Names PCPV\_gp124; Ankyrin/F-box protein

---

Official Symbol Ankyrin/F-box protein

---

Gene ID [8797325](#)

---

UniProt ID [D3IZR9](#)

---