

Tanapox, Yaba-like Disease Virus () 18L Recombinant Protein

Made to Order

Cat: MPYF-1122-KX136

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Tanapox, Yaba-like Disease Virus () 18L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	18L
Species	Tanapox, Yaba-like Disease Virus
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	18L
Full Name	18L protein
Alternative Names	18L; 18L protein; Ydvgp018
Gene ID	918730
UniProt ID	Q9DHU4