

Tanapox, Yaba-like Disease Virus () 12L Recombinant Protein

Made to Order

Cat: MPYF-1122-KX164

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Tanapox, Yaba-like Disease Virus () 12L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	12L
Species	Tanapox, Yaba-like Disease Virus
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	12L
Full Name	12L protein
Alternative Names	12L; 12L protein; Ydvgp012
Gene ID	918701
UniProt ID	Q9DHV0