

## Anti-Squirrel Parapoxvirus SQPV\_0150 Monoclonal Antibody (MP-K2507) (Mouse IgG)

Cat: MPYF-1122-KX1707

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Squirrel Parapoxvirus SQPV\_0150. The mAb with highest affinity will be selected. The selected antibody recognizes Squirrel Parapoxvirus SQPV\_0150. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

### Product Description

Target	SQPV_0150
Species Reactivity	Squirrel Parapoxvirus
Strain	Red squirrel UK
Cross Reactivity	SQPV_0150
Specificity	This antibody recognizes monkeypox virus (MPXV) SQPV_0150.
Antibody Isotype	Mouse IgG
Clone	MP-K2507
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

### Target

Target	SQPV_0150
Gene Name	SQPV_0150
Introduction	SQPV is a virus that causes the deadly disease squirrelpox in the United Kingdom red squirrels. The virus is commonly carried by gray squirrels, which rarely die from the disease. In red squirrels, the virus causes skin ulcers, lesions, and scabs. It is also accompanied by swelling near the eyes, mouth, feet, and genitals.
Alternative Names	SQPV_0150; hypothetical pox protein

Official Symbol            hypothetical pox protein

---

Gene ID                    [18158373](#)

---

UniProt ID                [Q1HTT7](#)

---