

Tanapox, Yaba-like Disease Virus () 114R Recombinant Protein

Made to Order

Cat: MPYF-1122-KX222

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Tanapox, Yaba-like Disease Virus () 114R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	114R
Species	Tanapox, Yaba-like Disease Virus
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	114R
Full Name	114R protein
Alternative Names	114R; 114R protein; Ydvgp115
Gene ID	918641
UniProt ID	Q9DHJ9