

## Tanapox, Yaba-like Disease Virus ( ) 106L Recombinant Protein

Made to Order

Cat: MPYF-1122-KX244

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Tanapox, Yaba-like Disease Virus ( ) 106L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

Target Protein	106L
Species	Tanapox, Yaba-like Disease Virus
Protein Class	Virology

### Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### Target

Target Protein	106L
Full Name	106L protein
Alternative Names	106L; 106L protein; Ydvgp107
Gene ID	<a href="#">918618</a>
UniProt ID	<a href="#">Q9DHK7</a>