

Yaba Monkey Tumor Virus (VR587) YMTVg21L Recombinant Protein

Made to Order

Cat: MPYF-1122-KX289

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Yaba Monkey Tumor Virus (VR587) YMTVg21L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	YMTVg21L
Species	Yaba Monkey Tumor Virus
Strain	VR587
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	YMTVg21L
Full Name	putative transmembrane protein
Alternative Names	YMTVg21L; putative transmembrane protein; integral membrane proteins that spans the internal and external surface of the lipid bilayer; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	2943707

