

Sealpox Virus (AFK76s1) CGV03_gp073 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX32

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sealpox Virus (AFK76s1) CGV03_gp073 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CGV03_gp073
Species	Sealpox Virus
Strain	AFK76s1
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CGV03_gp073
Full Name	RNA polymerase subunit RPO19-like protein
Alternative Names	CGV03_gp073; RNA polymerase subunit RPO19-like protein; SePPVgORF073
Gene ID	33349914
UniProt ID	A0A1Z3GCR7