

Sealpox Virus (AFK76s1) CGV03_gp061 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX42

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sealpox Virus (AFK76s1) CGV03_gp061 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CGV03_gp061
Species	Sealpox Virus
Strain	AFK76s1
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CGV03_gp061
Full Name	putative RNA polymerase subunit RPO18-like protein
Alternative Names	CGV03_gp061; putative RNA polymerase subunit RPO18-like protein; SePPVgORF061
Gene ID	33349902
UniProt ID	A0A1Z3GCS8