

Sealpox Virus (AFK76s1) CGV03_gp060 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX43

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sealpox Virus (AFK76s1) CGV03_gp060 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CGV03_gp060
Species	Sealpox Virus
Strain	AFK76s1
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CGV03_gp060
Full Name	putative early transcription factor VETF-like protein
Alternative Names	CGV03_gp060; putative early transcription factor VETF-like protein; SePPVgORF060
Gene ID	33349901
UniProt ID	A0A1Z3GCU2