

Sheeppox Virus (TU-V02127) SPPV_97 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX455

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sheeppox Virus (TU-V02127) SPPV_97 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SPPV_97
Species	Sheeppox Virus
Strain	TU-V02127
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SPPV_97
Full Name	Virion core protein P4a
Alternative Names	SPPV_97; Virion core protein P4a
Gene ID	944682
UniProt ID	A0A3F2YKQ8