

Sheeppox Virus (TU-V02127) SPPV_73 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX472

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sheeppox Virus (TU-V02127) SPPV_73 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SPPV_73
Species	Sheeppox Virus
Strain	TU-V02127
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SPPV_73
Full Name	DNA topoisomerase type I
Alternative Names	SPPV_73; DNA topoisomerase type I; ATP-independent enzyme that regulates the number of topological links between two DNA strands; relaxes positively and/or negatively supercoiled DNA
Gene ID	944656
UniProt ID	A0A3F2YKN1