

Sheeppox Virus (TU-V02127) SPPV_28 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX488

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sheeppox Virus (TU-V02127) SPPV_28 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SPPV_28
Species	Sheeppox Virus
Strain	TU-V02127
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SPPV_28
Full Name	poly(A) polymerase large subunit
Alternative Names	SPPV_28; poly(A) polymerase large subunit; catalyses template-independent extension of the 3'-end of a DNA or RNA; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	944633

