

## Sealpox Virus (AFK76s1) CGV03\_gp052 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX50

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sealpox Virus (AFK76s1) CGV03\_gp052 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |               |
|----------------|---------------|
| Target Protein | CGV03_gp052   |
| Species        | Sealpox Virus |
| Strain         | AFK76s1       |
| Protein Class  | Virology      |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                                          |
|-------------------|--------------------------------------------------------------------------|
| Target Protein    | CGV03_gp052                                                              |
| Full Name         | late transcription factor VLTF-4-like protein                            |
| Alternative Names | CGV03_gp052; late transcription factor VLTF-4-like protein; SePPVgORF052 |
| Gene ID           | <a href="#">33349893</a>                                                 |
| UniProt ID        | <a href="#">A0A1Z4CGG4</a>                                               |