

## Sheeppox Virus (TU-V02127) SPPV\_40 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX545

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sheeppox Virus (TU-V02127) SPPV\_40 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

Target Protein	SPPV_40
Species	Sheeppox Virus
Strain	TU-V02127
Protein Class	Virology

### Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### Target

Target Protein	SPPV_40
Full Name	hypothetical protein
Alternative Names	SPPV_40; hypothetical protein; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the structural proteins
Gene ID	<a href="#">944639</a>

