

Goatpox Virus (Pellor) GTPV_gp103 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX652

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Goatpox Virus (Pellor) GTPV_gp103 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	GTPV_gp103
Species	Goatpox Virus
Strain	Pellor
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	GTPV_gp103
Full Name	hypothetical protein
Alternative Names	GTPV_gp103; hypothetical protein
Gene ID	5228340
UniProt ID	A0A1B2LPQ7