

Goatpox Virus (Pellor) GTPV_gp114 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX699

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Goatpox Virus (Pellor) GTPV_gp114 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	GTPV_gp114
Species	Goatpox Virus
Strain	Pellor
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	GTPV_gp114
Full Name	hypothetical protein
Alternative Names	GTPV_gp114; hypothetical protein
Gene ID	5228262
UniProt ID	A0A1B2LPR9