

Lumpy Skin Disease Virus (Neethling 2490) LSDV050 Recombinant Protein Made to Order

Cat: MPYF-1122-KX726

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Lumpy Skin Disease Virus (Neethling 2490) LSDV050 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	LSDV050
Species	Lumpy Skin Disease Virus
Strain	Neethling 2490
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	LSDV050
Full Name	putative metalloprotease
Alternative Names	LSDV050; putative metalloprotease; metalloendopeptidase that likely plays a role in the maturation of viral proteins; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	921658

