

Sealpox Virus (AFK76s1) CGV03_gp114 Recombinant Protein

Made to Order

Cat: MPYF-1122-KX8

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Sealpox Virus (AFK76s1) CGV03_gp114 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CGV03_gp114
Species	Sealpox Virus
Strain	AFK76s1
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CGV03_gp114
Full Name	interleukin-10-like protein
Alternative Names	CGV03_gp114; interleukin-10-like protein; SePPVgORF114
Gene ID	33349955
UniProt ID	A0A1Z3GCP3