

Anti-Volepox virus VPXV-CA-129 Monoclonal Antibody (MP-K1732) (Mouse IgG)

Cat: MPYF-1122-KX932

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Volepox virus VPXV-CA-129. The mAb with highest affinity will be selected. The selected antibody recognizes Volepox virus VPXV-CA-129. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	VPXV-CA-129
Species Reactivity	Volepox virus
Strain	CA
Cross Reactivity	VPXV-CA-129
Specificity	This antibody recognizes monkeypox virus (MPXV) VPXV-CA-129.
Antibody Isotype	Mouse IgG
Clone	MP-K1732
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	VPXV-CA-129
Gene Name	VPXV-CA-129
Introduction	The Orthopoxvirus family includes the double-stranded DNA virus VPX. On the chorioallantoic membranes of chick embryos, VPX caused pinpoint, nonhemorrhagic pocks. There was a significant amount of cross-hybridization between the DNAs of VPX and other orthopoxviruses, however HindIII cleavage site maps revealed distinct ce

ntnal and terminal regions between VPX (222.8 kbp) and mapped DNAs of other orthopoxviruses.

Alternative Names	BI079_gp129; viral membrane formation protein; VPXV-CA-129
Official Symbol	viral membrane formation protein
Gene ID	29064076
UniProt ID	A0A1C9KCF2