
Anti-Sealpox Virus SePPVgORF086 Monoclonal Antibody (MP-K2632) (Mouse IgG)

Cat: MPYF-1222-KX1166

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Sealpox Virus SePPVgORF086. The mAb with highest affinity will be selected. The selected antibody recognizes Sealpox Virus SePPVgORF086. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	SePPVgORF086
Species Reactivity	Sealpox Virus
Strain	AFK76s1
Cross Reactivity	SePPVgORF086
Specificity	This antibody recognizes Sealpox Virus SePPVgORF086.
Antibody Isotype	Mouse IgG
Clone	MP-K2632
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	SePPVgORF086
Gene Name	SePPVgORF086
Introduction	Sealpox virus belongs to a separate species within the genus Parapoxvirus, which induces common infection in pinniped species such as harbor seals from North America and grey seals from Canada and Europe. The sealpox-induced disease is characterized by lesions of the skin, teats, and mucosa of the oral cavity in those seals.
Alternative Names	CGV03_gp086; hypothetical protein; SePPVgORF086

Official Symbol hypothetical protein

Gene ID [33349927](#)

UniProt ID [A0A1Z4CGF9](#)
