

Anti-Yaba Monkey Tumor Virus YMTVg78R Monoclonal Antibody (MP-K2882) (Mouse IgG)

Cat: MPYF-1222-KX1416

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Yaba Monkey Tumor Virus YMTVg78R. The mAb with highest affinity will be selected. The selected antibody recognizes Yaba Monkey Tumor Virus YMTVg78R. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	YMTVg78R
Species Reactivity	Yaba Monkey Tumor Virus
Strain	VR587
Cross Reactivity	YMTVg78R
Specificity	This antibody recognizes Yaba Monkey Tumor Virus YMTVg78R.
Antibody Isotype	Mouse IgG
Clone	MP-K2882
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	YMTVg78R
Gene Name	YMTVg78R
Introduction	The Yaba monkey tumor virus is a type of DNA poxvirus belonging to the genus Yatapoxvirus. It was first found in a colony of rhesus monkeys in Yaba, Lagos, Nigeria in 1956. While many poxviruses are host specific, the YMTV can not only infect and ca

use cutaneous nodules in primates, including rhesus monkeys, baboons, and cynomolgus monkeys, but also infect humans and cause mild and localized skin disease.

Alternative Names	YMTVg78R; 78R
Official Symbol	78R
Gene ID	2943628
UniProt ID	Q6TUT7