

Anti-Sheepox Virus SPPV_128 Monoclonal Antibody (MP-K3019) (Mouse IgG)

Cat: MPYF-1222-KX1553

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Sheepox Virus SPPV_128. The mAb with highest affinity will be selected. The selected antibody recognizes Sheepox Virus SPPV_128. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	SPPV_128
Species Reactivity	Sheepox Virus
Strain	TU-V02127
Cross Reactivity	SPPV_128
Specificity	This antibody recognizes Sheepox Virus SPPV_128.
Antibody Isotype	Mouse IgG
Clone	MP-K3019
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	SPPV_128
Gene Name	SPPV_128
Introduction	The sheepox virus is a poxvirus belonging to the genus Capripoxvirus. The natural hosts of SPV are domesticated and wild sheep. The younger and older populations of the sheep and lactating sheep are more susceptible to SPV. SPV infection cause s a devastating viral systemic disease that is defined by constant fever, widespread skin pox lesions, vesicles on non-wool skin, generalized papules or nodules, interna

l lesions of the respiratory and gastrointestinal mucosa, and ultimately leads to death.

Alternative Names	SPPV_128; VARV B22R homologue
Official Symbol	VARV B22R homologue
Gene ID	944640
UniProt ID	A0A3F2YKS6