

Myxoma Virus (Lausanne) m023R Recombinant Protein

Made to Order

Cat: MPYF-1222-KX180

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Myxoma Virus (Lausanne) m023R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	m023R
Species	Myxoma Virus
Strain	Lausanne
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	m023R
Full Name	m023R
Alternative Names	m023R; MYXV_gp027
Gene ID	932193
UniProt ID	Q77PA5