

Myxoma Virus (Lausanne) m087L Recombinant Protein

Made to Order

Cat: MPYF-1222-KX239

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Myxoma Virus (Lausanne) m087L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	m087L
Species	Myxoma Virus
Strain	Lausanne
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	m087L
Full Name	m87L
Alternative Names	m087L; MYXV_gp091
Gene ID	932119
UniProt ID	Q9Q8L3