

Myxoma Virus (Lausanne) m150R Recombinant Protein

Made to Order

Cat: MPYF-1222-KX309

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Myxoma Virus (Lausanne) m150R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	m150R
Species	Myxoma Virus
Strain	Lausanne
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	m150R
Full Name	m150R
Alternative Names	m150R; MYXV_gp153
Gene ID	932103
UniProt ID	Q9Q8F4