

Swinepox (17077-99) SPV042 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX39

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Swinepox (17077-99) SPV042 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SPV042
Species	Swinepox
Strain	17077-99
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SPV042
Full Name	SPV042 ribonucleotide reductase large chain
Alternative Names	SPV042; SWPVgp042; SPV042 ribonucleotide reductase large chain
Gene ID	932440
UniProt ID	Q8V3Q2