

Shope Fibroma Virus (Rabbit Fibroma) (Kasza) SFV_s037L Recombinant Protein Made to Order

Cat: MPYF-1222-KX444

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Shope Fibroma Virus (Rabbit Fibroma) (Kasza) SFV_s037L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SFV_s037L
Species	Shope Fibroma Virus (Rabbit Fibroma)
Strain	Kasza
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SFV_s037L
Full Name	gp037L
Alternative Names	SFV_s037L; gp037L
Gene ID	1486880
UniProt ID	Q9Q934