

Shope Fibroma Virus (Rabbit Fibroma) (Kasza) SFV_s032R Recombinant Protein Made to Order

Cat: MPYF-1222-KX479

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Shope Fibroma Virus (Rabbit Fibroma) (Kasza) SFV_s032R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	SFV_s032R
Species	Shope Fibroma Virus (Rabbit Fibroma)
Strain	Kasza
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	SFV_s032R
Full Name	Hypothetical protein
Alternative Names	SFV_s032R; hypothetical protein; similar to vaccinia virus strain Copenhagen E6R; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a biphasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	1486875

