

Canarypox Virus (VR-111) CNPV007 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX515

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV007 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CNPV007
Species	Canarypox Virus
Strain	VR-111
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CNPV007
Full Name	CNPV007 ankyrin repeat protein
Alternative Names	CNPV007; CNPV007 ankyrin repeat protein
Gene ID	2700362
UniProt ID	Q6VZZ0