

Canarypox Virus (VR-111) CNPV085 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX543

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV085 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CNPV085
Species	Canarypox Virus
Strain	VR-111
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CNPV085
Full Name	CNPV085 putative RNA phosphatase
Alternative Names	CNPV085; CNPV085 putative RNA phosphatase
Gene ID	2700311
UniProt ID	Q6VZR2