

Canarypox Virus (VR-111) CNPV125 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX589

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV125 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CNPV125
Species	Canarypox Virus
Strain	VR-111
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CNPV125
Full Name	CNPV125 variola B22R-like protein
Alternative Names	CNPV125; CNPV125 variola B22R-like protein
Gene ID	2700227
UniProt ID	Q6VZM2