

## Canarypox Virus (VR-111) CNPV253 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX705

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV253 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                 |
|----------------|-----------------|
| Target Protein | CNPV253         |
| Species        | Canarypox Virus |
| Strain         | VR-111          |
| Protein Class  | Virology        |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                         |
|-------------------|---------------------------------------------------------|
| Target Protein    | CNPV253                                                 |
| Full Name         | CNPV253 putative IMV membrane virulence factor          |
| Alternative Names | CNPV253; CNPV253 putative IMV membrane virulence factor |
| Gene ID           | <a href="#">2700316</a>                                 |
| UniProt ID        | <a href="#">Q6VZ94</a>                                  |