

Canarypox Virus (VR-111) CNPV308 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX759

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV308 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CNPV308
Species	Canarypox Virus
Strain	VR-111
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CNPV308
Full Name	CNPV308 ankyrin repeat protein
Alternative Names	CNPV308; CNPV308 ankyrin repeat protein
Gene ID	2700167
UniProt ID	Q6VZ39