

Canarypox Virus (VR-111) CNPV231 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX765

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Canarypox Virus (VR-111) CNPV231 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CNPV231
Species	Canarypox Virus
Strain	VR-111
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CNPV231
Full Name	CNPV231 MyD116-like domain protein
Alternative Names	CNPV231; CNPV231 MyD116-like domain protein
Gene ID	2700149
UniProt ID	Q6VZB6