

Fowlpox Virus (NVSL) FPV157 Recombinant Protein

Made to Order

Cat: MPYF-1222-KX830

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Fowlpox Virus (NVSL) FPV157 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	FPV157
Species	Fowlpox Virus
Strain	NVSL
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	FPV157
Full Name	p28-like protein
Alternative Names	FPV157; p28-like protein; similar to Homo sapiens GB:U41315; contains a KIL-A-N domain (conserved DNA-binding domain) that is found in a wide range of proteins of large bacterial and eukaryotic DNA viruses; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins; canarypox and fowlpox are members of the Avipoxvirus genus

and infect nonmammalian hosts; Fowlpox virus infects chickens and turkeys while Canarypox has a broader host-range and infects passeriform (song) birds.

Gene ID	1486705
---------	-------------------------

UniProt ID	Q9J573
------------	------------------------
